

HIRU Start Guide

Updating system

Hiru | The first multi-platform
eye tracker in the world.



Updates

We have released a new version of Hiru.

For using Hiru with EasyClick, Grid 3 or iOS apps, Hiru has to be **updated in a Windows device**. Once updated, it can be used either in Windows or iOS.


For updating Hiru to the latest firmware, please **download the last version of HiruSystray** app.

In the following pages, an explanation of how to update Hiru is included and in the last page, the links for downloading our apps.


Installation of updates

Install HiruSystray from this link:

<https://downloads.irisbond.com/systray>




Click Next




Accept the terms and Install

Updates

Installation of updates




Click Finish




Search the application in your PC and select it.

Updates

Installation of updates



The HiruSystray icon will be opened in the Windows Systray panel.






If the Hiru is connected to the PC, the Systray will periodically check if there are firmware updates available. It is also possible to check for the updates by clicking in *Check Updates*

Updates


Installation of updates

If an update is available, a notification will show up. Please, **DO NOT remove the Hiru during this process.**

These are the notifications that you will receive.



If the HIRU is already updated, you will receive the following notification



Now the Hiru is ready for using it with the last version of EasyClick, Grid 3 and our iOS app.

Updates

Installation of EasyClick and Grid 3 installers



EasyClick: Our computer access software. EasyClick is now compatible with Duo and Hiru in the same app. Download this new version of EasyClick [here](#).



Microsoft HID: It is possible to use Windows Eye Control with Hiru, and that can be done using our HiruSystray app. Download HiruSystray in this [link](#).



iOS Hiru app: Download the new version of the app directly from the App Store. Find it in [this link](#).

For more information about Hiru and the usage of applications, please look at our [Quick Start Guide](#).