

# HIRU Start Guide

## Windows Eye Control - Microsoft HID


*Hiru* | The first multi-platform  
eye tracker in the world.




# Windows Eye Control

Install HiruSystray from this link:

<https://downloads.irisbond.com/systray>




Click Next




Accept the terms and Install

# Windows Eye Control

## Install HiruSystray




Click Finish



Search the application in your PC and select it.

# Windows Eye Control

## Install HiruSystray



The HiruSystray icon will be opened in the Windows Systray panel.

## Open Windows Eye Control

### Eye control

Eye control lets you use eye tracking technology to control the mouse, type using an on screen keyboard, and communicate with people using text-to-speech.

### Use eye control (preview)

Turn on eye control (EN-US keyboard layout only)



[Get more info about eye control](#)


### Privacy information

We don't save or share eye tracking data. To learn more about how we use your data, see the privacy statement.

Go to Settings → Ease of Access → Eye Control →  
Turn on the toggle


# Windows Eye Control

Open Windows Eye Control




Windows Eye Control Bar will open

Activate HID integration for start using HIRU with Windows Eye Control





Go to HiruSystray Start HID




# Windows Eye Control

## Calibrate HIRU



After calibrating, start controlling the Windows Eye Control bar with HIRU

Stop HID connection if using other apps integrated with HIRU





IRISBOND  
see it possible

[support@irisbond.com](mailto:support@irisbond.com)